



## Security Bulletin - Digital Output Security Bulletin

As of 31-December-2008

Contact Carestream Health, Inc. Service to make arrangements for your equipment

Color Key Code: **Unacceptable risk**

**Low risk, fix will be in next major release or rollup patch**

N/A or N/A due to hardening, but may be applicable to other products

Note: any patches not on this report are deemed not applicable to Carestream products.

Product Susceptibility Analysis			DV8900			DV8150		CMI1000		DV6800	DV5800	MIM
Vulnerability ID	Vulnerability Description	Date Analyzed	1 / 2	3A / 3B / 3C	3D	1.0.1	1.0.3	1.1	1.3	1.0	1.0	6.1.1 / 6.1 / 6.0
MS08-078	Security Update for Internet Explorer (960714)	12/18/2008	—	—	—	—	—	—	—	—	—	—
MS08-073	Cumulative Security Update for Internet Explorer (958215)	12/12/2008	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening
MS08-070	Vulnerabilities in Visual Basic 6.0 Runtime Extended Files (ActiveX Controls) Could Allow Remote Code	12/9/2008	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening	N/A due to Hardening
MS08-067	Vulnerability in Server Service Could Allow Remote Code Execution (958644)	10/24/2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008	Low Risk/HZ24Oct2008
MSSA956391	Cumulative Security Update of ActiveX Kill Bits	10/15/2008	N/A	N/A	N/A	N/A	N/A	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening



Product Susceptibility Analysis			DV8900			DV8150		CMI1000		DV6800	DV5800	MIM
Vulnerability ID	Vulnerability Description	Date Analyzed	1 / 2	3A / 3B / 3C	3D	1.0.1	1.0.3	1.1	1.3	1.0	1.0	6.1.1 / 6.1 / 6.0
MS08-047	Vulnerability in IPsec Policy Processing Could Allow Information Disclosure (953733)	8/13/2008	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
MS08-045	Cumulative Security Update for Internet Explorer (953838)	8/13/2008	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening	N/A due to hardening
MS08-037	Vulnerabilities in DNS Could Allow Spoofing (953230)	7/10/2008	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A DNS support with V2.0, late Q4 2008	N/A